# REP Band™

## Resistive Exercise Products

Instruction Guide • Safety Precautions • Practical Tips  
Manuel d'instruction • Précautions de sécurité • Conseils pratiques  
Guía de instrucciones • Precauciones de seguridad • Consejos prácticos  
Benutzeranweisungen • Vorsichtsgebote • Praktische Hinweise  
Guida istruttiva • Precauzioni per la sicurezza • Consigli pratici  
Bruksanvisning • Säkerhetsbeaktanden • Praktiska tips

No Powder. No Latex. No Odor.  
Pas de poudre. Pas de latex. Pas d'odeur.  
Sin polvo. Sin látex. Sin olor.  
Pulverfrei. Latexfrei. Geruchlos.  
Niente talco. Niente lattice. Niente odore.  
Inget pulver. Inget latex. Ingen lukt.

REP Bands are patent pending. • Les REP Bands sont sous brevet en instance. • Las bandas REP tienen patente pendiente. • Patentanmeldung für REP-Bänder läuft • Il brevetto delle cinghie REP è in attesa di registrazione. • REP-bandet är patentsökt.

Magister Corporation **M**

(423) 265-3574 • Fax: (423) 265-4581  
E-mail: magistercorp@compuserve.com

In Europe, contact; • En Europe, contactez: • En Europa, comuníquese con:  
In Europa wenden Sie sich an: • In Europa, contattate: • Adress i Europa:  
Magister (UK) Ltd. • Three Chimneys House  
Upper Farringdon, Alton • Hampshire GU34 3ED • England

REP Band™ Resistive Exercise Products — Bands, Tubing and Loops — are designed to provide safe, easy and convenient exercise for rehabilitation, fitness conditioning and training.

- ◆ REP Band products eliminate the messy powder that gets on everything and everybody and can irritate sensitive skin.
- ◆ REP Band products are made from a synthetic elastomer and CONTAIN NO LATEX. That's important for any one who has latex allergies or exhibits latex allergic reactions.
- ◆ REP Band products are odor-free. Gone is the objectionable rubber smell. Odor-free is a real advantage when working with children and older adults.

**CAUTION:** The instructions contained in this guide are for reference only and not intended to act as a prescription for exercise. Before beginning any exercise program, it is strongly recommended that you consult with your health care provider. REP Band products are not toys; children should be supervised at all times.

### Safety Precautions

If used improperly, REP Band products can cause injury. Never use REP Band products in a way that causes the band to snap toward the head and eyes. Prior to use, always examine the REP Band products for any flaws or damage. Do not use REP Band products with holes or small nicks which may cause the band to tear or break during exercise. Replace a damaged product with a new one. Keep REP Band products away from sharp objects to avoid accidental punctures and tears.

### Practical Tips

For exercises which require tying REP Band exercise bands into loops, use a half bow or knot. Always make sure that the knot or half bow is secure prior to beginning the exercise. To eliminate the need for knots, use REP Band LOOP.

Long socks may be worn during use to prevent the band from pulling at leg hair.

Do not lock the joints when exercising. Keep elbows and knees slightly bent at all times.

Breathe evenly while exercising — don't hold your breath!

Maintain proper body alignment while performing standing exercises. Keep your shoulders back and stand up straight.

### Storage And Care Of Bands

Store your REP Band products in a cool (21°C/70°F), dark area. When using bands in chlorine pools, rinse bands after each use and allow to air dry.

© 1998, Magister Corporation. All rights reserved.  
© 1998, Magister Corporation. Tous droits réservés.  
© 1998, Magister Corporation. Reservados todos los derechos.  
© 1998, Magister Corporation. Alle Rechte vorbehalten.  
© 1998, Magister Corporation. Tutti i diritti sono riservati.  
© 1998 Magister Corporation. Med ensamrätt.

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
- ◆ Stand on band with both feet, grasp ends of each band.
- ◆ Slightly bend knees.
- ◆ Keep back straight and stomach muscles tightened.
- ◆ Pull with each arm to height of shoulders, pointing elbows out and away from body as you do so.
- ◆ Hold and return to starting position.

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- ◆ Tie band in a loop around the leg of a chair or other stationary object.
- ◆ Stand up straight with shoulders back.
- ◆ Grasp the back of the chair with one hand and extend the other arm out to the side for balance.
- ◆ Press leg backward, maintaining slight knee flexion. (Do not lock the knee.)
- ◆ Hold this position, then return to starting position.

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- ◆ Sit on floor in a long sitting position.
- ◆ Hold both ends of the band in one hand.
- ◆ Hook band under one foot. Bend that knee to 90° or less for a starting position.
- ◆ Take up the slack of the band in your hand.
- ◆ Press leg forward, keeping knee slightly bent.
- ◆ Do not lock the knee.
- ◆ Hold this position, then return to starting position.

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- ◆ Take the long sitting position, sitting up straight.
- ◆ Grip one end of the band in each hand.
- ◆ Loop the band under the balls of both feet.
- ◆ Pull both ends of the band toward the torso, leading with the elbows.
- ◆ Keep the elbows close to the torso.
- ◆ Hold this position, then return to starting position.

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- ◆ Sit with forearm resting on thigh with palm facing up.
- ◆ Grasp one end of the band in the hand, anchor the other end under the foot.
- ◆ Raise hand up as far as possible.
- ◆ Hold this position, then return to starting position.

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- ◆ Stand in a split-stance with knees bent.
- ◆ Grab one end of the band in each hand.
- ◆ Start with the ends of the band just above the knee joint.
- ◆ Curl, keeping wrists straight.
- ◆ Hold this position, then return to starting position.

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- ◆ Tie band in loop.
- ◆ Standing on band with one foot, place other foot inside loop.
- ◆ Balance yourself in a slight forward lean against chair back.
- ◆ Curl leg inside the loop to bring foot off floor and bend knee to a 70° angle.
- ◆ Hold this position, then return to starting position.

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- ◆ Sit in a chair with one arm down at your side, elbow at 90°.
- ◆ Loop band around palm.
- ◆ Rotate body away from side where band is anchored.
- ◆ Keep upper body straight.
- ◆ Hold this position, then return to starting position.

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